


Spring clamps for laboratory flasks

Accessories for **SONOREX** ultrasonic baths in the laboratory



- neither floating nor canting of flasks
- fast and easy fixing to baskets and holders with mesh sizes up to 12.5 × 12.5 mm
- material: stainless steel



	EK 10	EK 25	EK 50	EK 100	EK 250
For laboratory flask	10 ml	25 ml	50 ml	100 ml	250 ml
Up to diameter	31 mm	42 mm	52 mm	65 mm	85 mm
Minimum diameter	23 mm	30 mm	35 mm	40 mm	55 mm
Code No.	051	053	055	057	3259